

TECHNICAL DATASHEET Polypectomy Snare
--

Type: Polypectomy Snare (disposable)

Article-No.*: * see Table below

Material:

Component	Appearance	Additives	Material	Applicable Standard
Handle	Solid	N/A	ABS	ISO 8283-4-1992
Drive wire	Wire	N/A	304	-
Out Tube	Tube	N/A	PTFE	-
Loop	Wire	N/A	304	-
Conducting Plug	Solid	N/A	Hpd59-1	-
Conductive Block	Solid	N/A	Hpd59-1	-
Reinforcing Tube	Solid	N/A	304	-
Boosting tube	Solid	N/A	304	-

Measurements*: * see Table below

Intended Use: The snare is indicated for use in the ligation, electro-cautery, resection and removal of tissues and polyps.

Supplied / characteristic:

Sterile



Single-Use



Quality Assurance:

Continuous control during production, final control and release.
 DIN EN ISO 13485
 Traceability by LOT numbers.

Packaging:

1 piece per package

Storage conditions:

Store in a dry, cool and clean place and away from direct sunlight, heat and/or chemical fumes.

Disposal: Normal disposal within the scope of the waste ordinances of the respective state.

***Measurements / Article-No.:**

(*Table)

Art. No.	Outer dia.(mm)	Total Length(cm)	Open Size (mm)	Type
BP-40201	1.8	100	10	Oval
BP-40202	1.8	100	20	Oval
BP-40203	1.8	160	10	Oval
BP-40204	1.8	160	20	Oval
BP-40205	2.3	160	15	Oval
BP-40206	2.3	160	20	Oval
BP-40207	2.3	160	25	Oval
BP-40208	2.3	160	30	Oval
BP-40209	2.3	230	15	Oval
BP-40210	2.3	230	20	Oval
BP-40211	2.3	230	25	Oval
BP-40212	2.3	230	30	Oval
Art. No.	Outer dia.(mm)	Total Length(cm)	Open Size (mm)	Type
BP-40213	1.8	100	10	Hexagonal
BP-40214	1.8	100	20	Hexagonal
BP-40215	1.8	160	10	Hexagonal
BP-40216	1.8	160	20	Hexagonal
BP-40217	2.3	160	15	Hexagonal
BP-40218	2.3	160	20	Hexagonal
BP-40219	2.3	160	25	Hexagonal
BP-40220	2.3	160	30	Hexagonal
BP-40221	2.3	230	15	Hexagonal
BP-40222	2.3	230	20	Hexagonal
BP-40223	2.3	230	25	Hexagonal
BP-40224	2.3	230	30	Hexagonal
Art. No.	Outer dia.(mm)	Total Length(cm)	Open Size (mm)	Type
BP-40225	1.8	100	10	Crescent
BP-40226	1.8	100	20	Crescent
BP-40227	1.8	160	10	Crescent
BP-40228	1.8	160	20	Crescent
BP-40229	2.3	160	15	Crescent
BP-40230	2.3	160	20	Crescent
BP-40231	2.3	160	25	Crescent
BP-40232	2.3	160	30	Crescent
BP-40233	2.3	230	15	Crescent
BP-40234	2.3	230	20	Crescent
BP-40235	2.3	230	25	Crescent
BP-40236	2.3	230	30	Crescent