

TECHNICAL DATASHEET Cytology Brush

Type: Cytology Brush (disposable)

Article-No.*: * see Table below

Material:

Component	Appearance	Additives	Generic Name	Criterion for Implementation
Handle	Solid	N/A	ABS	ISO 8283-4-1992
Drive wire	Wire	N/A	304	-
Out Tube	Tube	N/A	PTFE	-
Bristle	Bristle	N/A	Nylon	-

Measurements *: * see Table below

Intended Use: The Cytology Brushes are used to collect cells from the bronchi and upper and lower gastrointestinal tracts. This device is supplied sterile and intended for single use only.

Supplied /characteristic:

Sterile



Single-Use



Quality Assurance: Continuous control during production, final control and release.
DIN EN ISO 13485
Traceability by LOT numbers.

Packaging: 1 piece per package

Storage conditions: Store in a dry, cool and clean place and away from direct sunlight, heat and/or chemical fumes.

Disposal: Normal disposal within the scope of the waste ordinances of the respective state.

***Measurements / Article-No.:**

(*Table)

Art. No.	Tube dia.(mm)	Total Length(cm)	Brush Type
BP-40801	1.8	105	Straight
BP-40802	1.8	120	Straight
BP-40803	2.3	160	Straight
BP-40804	2.3	180	Straight
BP-40805	2.3	230	Straight
Art. No.	Tube dia.(mm)	Total Length(cm)	Brush Type
BP-40806	1.8	105	U shape
BP-40807	1.8	120	U shape
BP-40808	2.3	160	U shape
BP-40809	2.3	180	U shape
BP-40810	2.3	230	U shape