

<b>TECHNICAL DATASHEET</b> <b>Biopsy Valve</b>
---

**Type:** Biopsy Valve (Disposable)

**Article-No.\*:** \* see Table below

Article-No.	Color	Channel
BP-40601	Black	For Olympus and Fujinon Endoscope use
BP-40602	Black	For Pentax Endoscope use

**Material:**

Component	Appearance	Additives	Generic Name	Criterion for Implementation
Biopsy Valve	Solid	N/A	Silicon dioxide	-

**Intended Use:** Disposable Biopsy valves provide access for endoscopic device passage and exchange, helps maintain insufflation and minimizes leakage of biomaterial from the biopsy port throughout the gastrointestinal endoscopic procedure.

**Supplied /characteristic:**

non-Sterile



Single-Use



**Quality Assurance:** Continuous control during production, final control and release.  
DIN EN ISO 13485  
Traceability by LOT numbers.

**Packaging:** 1 piece per package

**Storage conditions:** Store in a dry, cool and clean place and away from direct sunlight, heat and/or chemical fumes.

**Disposal:** Normal disposal within the scope of the waste ordinances of the respective state.