


Speciality Long Luer Needles

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION



Standard luer connector

Colour coded for easy identification of size

Colour coded dot corresponds to bevel

Polycarbonate hub prevents leakage

17° Quincke bevel

Speciality long luer needles compliant with ISO80369-7

Range of sizes and lengths

LUER

Buckingham

BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000
Fax: 023 8066 9002
info@blueboxmedical.co.uk
www.blueboxmedical.co.uk



PS060

04/17 - A

SPECIALITY NEEDLES

SPECIALITY LONG LUER NEEDLES

Reference	Colour Code	Gauge	Length	Units/Box
HTS1401	Pink	18G	50mm	25
HTS1404	Pink	18G	90mm	25
HTS1405	Pink	18G	125mm	25
HTS1601	Yellow	20G	60mm	25
HTS1603	Yellow	20G	90mm	25
HTS1606	Yellow	20G	125mm	25
HTS1607	Yellow	20G	150mm	25
HTS1608	Yellow	20G	180mm	25
HTS1801	Black	22G	38mm	25
HTS1803	Black	22G	90mm	25
HTS1807	Black	22G	150mm	25
HTS1809	Black	22G	180mm	25
HTS2003	Black	22G	50mm	25
HTS2006	Black	22G	120mm	25
HTS2005	Orange	25G	90mm	25