


Medi-Inn[®] Biopsy & Grasping Forceps

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION

ALLIGATOR OR OVAL CUPS

Available coated and non-coated, with or without spike



High precision cut avoids risk of tissue damage

Integrated ergonomic handle for comfortable use

Includes tissue sampler

Length marking on coated version


GRASPING FORCEPS WITH PRONGS

For removal of small and large polyps

Large opening width

Atraumatic distal endions reduce tissue damage to a minimum

Integrated ergonomic handle for comfortable use



3, 4 and 5 prong



BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002
info@blueboxmedical.co.uk
www.blueboxmedical.co.uk



PS062

10/17 - B

ENDOSCOPY

MEDI-INN® BIOPSY & GRASPING FORCEPS

Reference (without spike)	NPC	Reference (with spike)	NPC	Outer Diameter	Working Length	Case Qty
BIOPSY FORCEPS (NON-COATED - OVAL CUPS)						
BP-40101	FFI1535	BP-40108	FFI1542	1.8mm	100cm	10 x 10
BP-40102	FFI1536	BP-40109	FFI1543	1.8mm	120cm	10 x 10
BP-40103	FFI1537	BP-40110	FFI1560	1.8mm	150cm	10 x 10
BP-40104	FFI1538	BP-40111	FFI1561	2.3mm	160cm	10 x 10
BP-40105	FFI1539	BP-40112	FFI1562	2.3mm	180cm	10 x 10
BP-40106	FFI1540			2.3mm	210cm	10 x 10
BP-40107	FFI1541	BP-40114	FFI1544	2.3mm	230cm	10 x 10
BIOPSY FORCEPS (NON-COATED - ALLIGATOR CUPS)						
BP-40129	FFI1559	BP-40136	FFI1572	1.8mm	100cm	10 x 10
BP-40130	FFI1567	BP-40137	FFI1573	1.8mm	120cm	10 x 10
BP-40131	FFI1568	BP-40138	FFI1574	1.8mm	150cm	10 x 10
-	-	BP-40139	FFI1575	2.3mm	160cm	10 x 10
BP-40133	FFI1569	BP-40140	FFI1576	2.3mm	180cm	10 x 10
BP-40134	FFI1570	BP-40141	FFI1577	2.3mm	210cm	10 x 10
BP-40135	FFI1571	BP-40142	FFI1578	2.3mm	230cm	10 x 10
BIOPSY FORCEPS (COATED - OVAL CUPS)						
BP-40115	FFI1545	BP-40122	FFI1552	1.8mm	100cm	10 x 10
BP-40116	FFI1546	BP-40123	FFI1553	1.8mm	120cm	10 x 10
BP-40117	FFI1547	BP-40124	FFI1554	1.8mm	150cm	10 x 10
BP-40118	FFI1548	BP-40125	FFI1555	2.3mm	160cm	10 x 10
BP-40119	FFI1549	BP-40126	FFI1556	2.3mm	180cm	10 x 10
BP-40120	FFI1550	BP-40127	FFI1557	2.3mm	210cm	10 x 10
BP-40121	FFI1551	BP-40128	FFI1558	2.3mm	230cm	10 x 10
BIOPSY FORCEPS (COATED - ALLIGATOR CUPS)						
BP-40143	FFI1579	-	-	1.8mm	100cm	10 x 10
BP-40144	FFI1580	BP-40151	FFI1586	1.8mm	120cm	10 x 10
BP-40145	FFI1581	BP-40152	FFI1587	1.8mm	150cm	10 x 10
BP-40146	FFI1582	BP-40153	FFI1588	2.3mm	160cm	10 x 10
BP-40147	FFI1583	BP-40154	FFI1589	2.3mm	180cm	10 x 10
BP-40148	FFI1584	BP-40155	FFI1590	2.3mm	210cm	10 x 10
BP-40149	FFI1585	BP-40156	FFI1591	2.3mm	230cm	10 x 10

Reference	NPC	Tube outer ø	Working length	Case Qty
3-PRONG GRASPING FORCEPS				
BP-40909	FFI1564	2.3mm	160cm	10 x 10
BP-40912	FFI1594	2.3mm	230cm	10 x 10
4-PRONG GRASPING FORCEPS				
BP-40910	FFI1592	2.3mm	160cm	10 x 10
BP-40913	FFI1563	2.3mm	230cm	10 x 10
5-PRONG GRASPING FORCEPS				
BP-40911	FFI1593	2.3mm	160cm	10 x 10
BP-40914	FFI1566	2.3mm	230mm	10 x 10



BLUE BOX
MEDICAL LTD

Tel: 023 8066 9000 Fax: 023 8066 9002
info@blueboxmedical.co.uk
www.blueboxmedical.co.uk