


Medi-Inn[®] Polypectomy Snares

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION



ENDOSCOPY



BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002
info@blueboxmedical.co.uk
www.blueboxmedical.co.uk



PS063

10/17 - B

MEDI-INN® POLYPECTOMY SNARES

Reference	NPC	Shape	Tube Outer ø	Working Length	Opening	Case Qty
BP-40201	FBE22783	Oval	1.8mm	100cm	10mm	10 x 100
BP-40202	FBE22794	Oval	1.8mm	100cm	20mm	10 x 100
BP-40203	FBE22795	Oval	1.8mm	160cm	10mm	10 x 100
BP-40204	FBE22796	Oval	1.8mm	160cm	20mm	10 x 100
BP-40205	FBE22797	Oval	2.3mm	160cm	15mm	10 x 100
BP-40206	FBE22798	Oval	2.3mm	160cm	20mm	10 x 100
BP-40207	FBE22799	Oval	2.3mm	160cm	25mm	10 x 100
BP-40208	FBE22800	Oval	2.3mm	160cm	30mm	10 x 100
BP-40209	FBE22801	Oval	2.3mm	230cm	15mm	10 x 100
BP-40210	FBE22802	Oval	2.3mm	230cm	20mm	10 x 100
BP-40211	FBE22803	Oval	2.3mm	230cm	25mm	10 x 100
BP-40212	FBE22804	Oval	2.3mm	230cm	30mm	10 x 100
BP-40213	FBE22805	Hexagonal	1.8mm	100cm	10mm	10 x 100
BP-40214	FBE22806	Hexagonal	1.8mm	100cm	20mm	10 x 100
BP-40215	FBE22807	Hexagonal	1.8mm	160cm	10mm	10 x 100
BP-40216	FBE22808	Hexagonal	1.8mm	160cm	20mm	10 x 100
BP-40217	FBE22809	Hexagonal	2.3mm	160cm	15mm	10 x 100
BP-40218	FBE22810	Hexagonal	2.3mm	160cm	20mm	10 x 100
BP-40219	FBE22811	Hexagonal	2.3mm	160cm	25mm	10 x 100
BP-40220	FBE22812	Hexagonal	2.3mm	160cm	30mm	10 x 100
BP-40221	FBE22813	Hexagonal	2.3mm	230cm	15mm	10 x 100
BP-40222	FBE22814	Hexagonal	2.3mm	230cm	20mm	10 x 100
BP-40223	FBE22815	Hexagonal	2.3mm	230cm	25mm	10 x 100
BP-40224	FBE22816	Hexagonal	2.3mm	230cm	30mm	10 x 100
BP-40225	FBE22817	Crescent	1.8mm	100cm	10mm	10 x 100
BP-40226	FBE22818	Crescent	1.8mm	100cm	20mm	10 x 100
BP-40227	FBE22819	Crescent	1.8mm	160cm	10mm	10 x 100
BP-40228	FBE22820	Crescent	1.8mm	160cm	20mm	10 x 100
BP-40229	FBE22821	Crescent	2.3mm	160cm	15mm	10 x 100
BP-40230	FBE22822	Crescent	2.3mm	160cm	20mm	10 x 100
BP-40231	FBE22823	Crescent	2.3mm	160cm	25mm	10 x 100
BP-40232	FBE22824	Crescent	2.3mm	160cm	30mm	10 x 100
BP-40233	FBE22825	Crescent	2.3mm	230cm	15mm	10 x 100
BP-40234	FBE22826	Crescent	2.3mm	230cm	20mm	10 x 100
BP-40235	FBE22827	Crescent	2.3mm	230cm	25mm	10 x 100



BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002

info@blueboxmedical.co.uk

www@blueboxmedical.co.uk