

# Medi-Inn<sup>®</sup> Cytology Brushes

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION


Metal ball head for optimum x-ray visibility

Brush smear for cytological analysis

Microbores for better diagnostic results

Protection of the cytologic smear by proven tube closure

Integrated ergonomic handle for comfortable use



ENDOSCOPY



**BLUE BOX**  
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002  
info@blueboxmedical.co.uk  
www.blueboxmedical.co.uk



PS064

10/17 - B

## MEDI-INN® CYTOLOGY BRUSHES

Reference	NPC	Shape	Tube Outer ø	Working Length	Case Qty
BP-40801	IMC600	Straight	1.8mm	105cm	100
BP-40802	IMC000	Straight	1.8mm	120cm	100
BP-40803	IMC007	Straight	2.3mm	160cm	100
BP-40804	IMC002	Straight	2.3mm	180cm	100
BP-40806	IMC003	U-Shape	1.8mm	105cm	100
BP-40807	IMC001	U-Shape	1.8mm	120cm	100
BP-40808	IMC004	U-Shape	2.3mm	160cm	100
BP-40809	IMC005	U-Shape	2.3mm	180cm	100



**BLUE BOX**   
MEDICAL LTD

Tel: 023 8066 9000 Fax: 023 8066 9002  
info@blueboxmedical.co.uk  
www@blueboxmedical.co.uk