


# Medi-Inn<sup>®</sup> Foreign Body Retrieval Nets

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION

Removal of difficult foreign bodies

Excellent visibility and high precision manoeuvrability



Strong mesh material giving high security to collect foreign bodies as well as polyps and polyp fragments

Integrated ergonomic handle for comfortable use

ENDOSCOPY



BLUE BOX  
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002  
info@blueboxmedical.co.uk  
www.blueboxmedical.co.uk



PS065

10/17 - B

**MEDI-INN® FOREIGN BODY RETRIEVAL NETS**

Reference	NPC	Tube Outer ø	Working Length	Opening	Case Qty
BP-40901	FBE22792	1.8mm	160cm	20mm	10 x 10
BP-40902	FBE22787	1.8mm	160cm	25mm	10 x 10
BP-40903	FBE22782	1.8mm	160cm	30mm	10 x 10
BP-40904	FBE22793	1.8mm	160cm	40mm	10 x 10
BP-40905	FBE22789	2.5mm	230cm	20mm	10 x 10
BP-40906	FBE22790	2.5mm	230cm	25mm	10 x 10
BP-40907	FBE22791	2.5mm	230cm	30mm	10 x 10
BP-40908	FBE22786	2.5mm	230cm	40mm	10 x 10



**BLUE BOX**   
MEDICAL LTD

Tel: 023 8066 9000 Fax: 023 8066 9002  
info@blueboxmedical.co.uk  
www@blueboxmedical.co.uk