

Medi-Inn[®] Cleaning Brushes

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION

Ball heads on brush ends and rounded brush tips for gentle cleaning

Single/dual-ended cleaning brush made of soft, resistant nylon

Combination brush (dual-ended) for air/water valve and working channel

Brush head diameter allows maximum mechanical cleaning of the endoscope channel

Sweeper brush for advanced cleaning

Combination brush (dual-ended Xmas) for air/water valve and working channel

Conical brush head for optimised cleaning

Multipoint dual-ended cleaning brush combines the advantages of a dual-ended brush with a very efficient sweeper

Brush made of soft resistant nylon

For air and water channel

Rounded brush head for gentle cleaning

Vinyl cap on the end of the handle

Effective and gentle cleaning performance with U-shape



BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002
info@blueboxmedical.co.uk
www.blueboxmedical.co.uk



PS067

10/17 - B

ENDOSCOPY

MEDI-INN® CLEANING BRUSHES

Reference	NPC	Colour	Catheter		Brush 1		Brush 2		Case Qty
			∅	Length	∅	Length	∅	Length	
SINGLE-/DUAL-ENDED									
BP-41450	IMC008	Yellow	1.7mm	230cm	5.0mm	20mm	-	-	4 x 100
BP-41469	MJS507	Yellow	1.7mm	230cm	6.0mm	20mm	-	-	4 x 100
BP-41449	MJS540	Yellow	1.0mm	100cm	2.0mm	15mm	-	-	4 x 100
BP-40399	MJS506	Blue	1.0mm	130cm	1.5mm	10mm	-	-	4 x 100
BP-41454	MJS517	Yellow	1.7mm	180cm	5.0mm	20mm	5.0mm	20mm	4 x 100
BP-40392	MJS518	Yellow	1.7mm	230cm	5.0mm	20mm	5.0mm	20mm	4 x 100
BP-41455	MJS519	Blue	1.7mm	230cm	6.0mm	20mm	6.0mm	20mm	4 x 100
BP-41472	MJS520	Blue	1.7mm	300cm	6.0mm	20mm	6.0mm	20mm	4 x 100
BP-41456	MJS521	Yellow	1.8mm	230cm	6.0mm	20mm	6.0mm	20mm	4 x 100
BP-41470	MJS514	Green	1.5mm	100cm	3.0mm	20mm	3.0mm	20mm	4 x 100
BP-41471	MJS515	Green	1.5mm	180cm	3.0mm	20mm	3.0mm	20mm	4 x 100
BP-41453	MJS516	Green	1.5mm	230cm	3.0mm	20mm	3.0mm	20mm	4 x 100
BP-41473	MJS509	Yellow	1.0mm	120cm	2.0mm	15mm	2.0mm	15mm	4 x 100
BP-41451	MJS508	Blue	1.0mm	130cm	1.5mm	10mm	1.5mm	10mm	4 x 100
BP-41452	MJS513	Blue	1.0mm	230mm	2.0mm	15mm	2.0mm	15mm	4 x 100

Reference	NPC	Colour	Catheter		Valve Brush		Channel Brush		Case Qty
			∅	Length	∅	Length	∅	Length	
COMBINATION BRUSHES - DUAL-ENDED CHANNEL & VALVE									
BP-40398	MJS525	White	1.7mm	230cm	10mm	40mm	5mm	15mm	4 x 100
BP-40395	MJS524	Sky Blue	1.7mm	230cm	10mm	40mm	6mm	20mm	4 x 100
BP-41462	MJS523	Blue	1.7mm	230cm	11mm	40mm	6mm	20mm	4 x 100

Reference	NPC	Colour	Catheter		Brush		Sweeper		Case Qty
			∅	Length	∅	Length	∅	Length	
SWEEPER BRUSHES									
BP-40377	MJS526	Yellow	1.7mm	230cm	5mm	20mm	5mm	300mm	4 x 100
BP-40378	MJS527	Blue	1.7mm	160cm	5mm	20mm	5mm	300mm	4 x 100
BP-40379	MJS528	Orange	1.5mm	160cm	3mm	20mm	3mm	300mm	4 x 100
BP-40380	MJS510	Green	1.5mm	120cm	3mm	20mm	3mm	300mm	4 x 100

Reference	NPC	Colour	Catheter		Valve Brush		Sweeper		Channel		Case Qty
			∅	Length	∅	Length	∅	Length	∅	Length	
COMBINATION BRUSH - DUAL-ENDED (XMAS)											
BP-41468	MJS512	Sky Blue	1.7mm	230cm	7 - 4.5mm	20mm	5mm	100mm	5 - 3mm	10mm	4 x 100

Reference	NPC	Colour	Catheter		Brush 1		Sweeper		Brush 2		Case Qty
			∅	Length	∅	Length	∅	Length	∅	Length	
MULTIPOINT DUAL-ENDED CLEANING BRUSH											
BP-40393	MJS522	Blue	1.7mm	230cm	6mm	20mm	6mm	20mm	6mm	20mm	4 x 100

Reference	NPC	Colour	Length	Brush		Case Qty
				∅	Length	
U-SHAPE BRUSH						
BP-40335	MJS511	Blue	15cm	70mm	25mm	4 x 30

Reference	NPC	Colour	Length	Valve Brush 1		Valve Brush 2		Case Qty
				∅	Length	∅	Length	
SINGLE-/DUAL-ENDED VALVE BRUSHES								
BP-40387	FGY063	Yellow	15cm	11mm	35mm	5mm	20mm	4 x 200
BP-40388	FGY064	Sky Blue	15cm	11mm	35mm	5mm	20mm	4 x 200
BP-40389	FGY065	Dark Blue	15cm	11mm	35mm	5mm	20mm	4 x 200
BP-40390	FGY066	Red	15cm	11mm	35mm	-	-	4 x 200
BP-40391	FGY067	Green	15cm	11mm	35mm	-	-	4 x 200