

Medi-Inn[®] Accessories

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION

POLYP TRAPS



Numbered chambers to identify the collected polyps and allow exact examination

Safety sleeve prevents tissue slipping

Flexible connection tubes with adaptors for suction pump

DIGNITY SHORTS



Universal size

Elastic waist and slit back

Breathable, non-woven fabric

Easy to open and close

Safe protection of patient and nursing staff

Optimal seal



Secured safety on the biopsy channel opening

IRRIGATION - WATER BOTTLE CAPS, TUBING & CONNECTOR

24-hour use of water bottle caps for prevention of cross contamination; with CO2 connector option; for all existing NaCl bottles

Irrigation scope connector for use with flushing hose in conjunction with flushing pump and rinsing fluids

24 hour use on tubing; single-use in order to support patient safety



ENDOSCOPE DISTAL PROTECTION & CLEANING SPONGES



Sponges for cleaning and for protection of distal tip of endoscopes



BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002

info@blueboxmedical.co.uk

www.blueboxmedical.co.uk

ENDOSCOPY



PS068

10/17 - B

MEDI-INN® ACCESSORIES

Reference	NPC	Description	Case Qty
BP-40701	FBE22785	Polyp Trap 4 Chambers	10 x 10
BP-40702	FBE22788	Polyp Trap 5 Chambers	10 x 10

Reference	NPC	Description	Case Qty
BP-40601	FGY068	Biopsy Valves for Olympus and Fujinon Series 5	4 x 200
BP-40602	FGY069	Biopsy Valves for Pentax Endoscopes	4 x 200

Reference	NPC	Description	Case Qty
BP-41501	FAL17690	Dignity Shorts	100

Reference	NPC	Description	Case Qty
BP-41201	FAL17691	Water Bottle Caps for NaCl Bottles, with CO2 - for Olympus	4 x 10
BP-41202	FAL17692	Water Bottle Caps for NaCl Bottles, with CO2 - for Fujinon	4 x 10
BP-41203	FAL17693	Water Bottle Caps for NaCl Bottles, with CO2 - for Pentax	4 x 10
BP-41205	FBE22784	Irrigation Tube for Olympus Pumps	4 x 10

Reference	NPC	Colour	Description	Diameter	Length	Outer ø	Inner ø	Case Qty
BP-40501	MJS502	Blue	Sponge protector	> 5mm	140mm	37mm	9-24mm	4 x 50
BP-40502	MJS503	Red	Sponge protector	< 5mm	140mm	37mm	4-24mm	4 x 50
BP-40503	MJS504	White	3-in-1 sponge protector & cleaner	>/< 5mm	140mm	37mm	-	4 x 50
BP-40504	MJS505	Green	Cleaning sponge	>/< 5mm	140mm	37mm	-	4 x 50



BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002
info@blueboxmedical.co.uk
www@blueboxmedical.co.uk