

AutoFuser® Disposable Silicone Balloon Infuser

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION

Continuous infusion for antibiotic, chemotherapy and pain management therapies



Easy fill

Optional PCA bolus

Multiple reservoir sizes

Self-priming

Hard, protective shell

Latex free

Carrying pouch with strap for ambulatory flexibility

BLUE BOX 
MEDICAL LTD

Tel: 023 8066 9000 Fax: 023 8066 9002

info@blueboxmedical.co.uk

www.blueboxmedical.co.uk

ACEMEDICAL



PS021

01/17 - B

PATIENT CARE

AUTO FUSER® DISPOSABLE SILICONE BALLOON INFUSER

Reference	Sizes		Bolus	Lock-out Time	Units/Box
	Volume	Flow Rate			
AA1505	150ml	5ml/hr	-	-	10
AA1507	150ml	2.5ml/hr	-	-	10
AA1900	100ml	0.5ml/hr	-	-	10
AA1902	100ml	2ml/hr	-	-	10
AA1905	100ml	5ml/hr	-	-	10
AA2002	275ml	2ml/hr	-	-	10
AA2005	275ml	5ml/hr	-	-	10
AA2008	275ml	8ml/hr	-	-	10
AA2010	275ml	10ml/hr	-	-	10
AA2021	275ml	2ml/hr	1ml	15 mins	10
AA2029	275ml	5ml/hr	2ml	30 mins	10
AA4004	400ml	4ml/hr	-	-	10
AA4005	400ml	5ml/hr	-	-	10

This represents a small sample of our comprehensive range of volumes, flow rates and bolus dosage. Please contact us for further details.

BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002

info@blueboxmedical.co.uk

www@blueboxmedical.co.uk

ACEMEDICAL