

Micro-Tip Atraumatic Spinal Needles

PLEASE SEE OVERLEAF FOR PRODUCT INFORMATION

Siliconised, medical quality stainless steel reduces friction when inserting the needle

New 'dual diameter' design combines the stiffness of a large gauge needle with a small needle tip

Reduction of dura mater trauma and PDPH rate consistent with needle tip gauge (27G)

Quick flush-back of spinal fluid reduces technique failure

Good flexibility of the tip reduces the risk of bone injury; designed to bend 180° without breaking

Greater internal diameter along most of the needle length improves ease of injection

360° magnifying window in the hub allows quicker recognition of spinal fluid

Centimetre marking on the needle tubing improves perception of the depth of penetration

Double colour code on the stylet hub reminds the user of the double diameter of the needle



Atraumatic, **conical-elliptic** tip



BLUE BOX
MEDICAL LTD

Unit 29, New Forest Enterprise Centre,
Chapel Lane, Totton, Southampton SO40 9LA



Tel: 023 8066 9000 Fax: 023 8066 9002

info@blueboxmedical.co.uk

www.blueboxmedical.co.uk



PS035

11/14 -A

SPINAL/EPIDURAL ANAESTHESIA

MICRO-TIP ATRAUMATIC SPINAL NEEDLES

Reference	Colour Code		Sizes		Introducer	Units/Box
			Gauge	Length		
27922-27	Black	Grey	22/27G	70mm	25mm	25
29922-27	Black	Grey	22/27G	90mm	30mm	25
21922-27	Black	Grey	22/27G	103mm	38mm	25
31922-27	Black	Grey	22/27G	140mm	42mm	25



BLUE BOX
MEDICAL LTD

Unit 29, New Forest Enterprise Centre,
Chapel Lane, Totton, Southampton SO40 9LA



Tel: 023 8066 9000 Fax: 023 8066 9002

info@blueboxmedical.co.uk

www.blueboxmedical.co.uk