


Urodynamics Catheters

BLADDER PRESSURE CATHETER

Holes cut at 1cm
and 2cm from
catheter tip


4.5fg x 2m
blue tint tube



BLADDER FILLING CATHETER

Hole cut between 10mm
and 15mm from
catheter tip


6fg/8fg/10fg/12fg x 385mm
clear PVC tube



DOUBLE LUMEN CYSTOMETRY CATHETER

Holes cut at 5mm
and 15mm from
catheter tip



6fg/8fg/10fg x 70cm
mid blue double lumen tube



RECTAL BALLOON CATHETER

PVC balloon
(with slit)

4.5fg x 2m
red tint tube



CATHETERS

Reference	Designation	Units/Box
BCC-045	Bladder Pressure Catheter 4.5fg	25
BCF-060	Bladder Filling Catheter 6fg	25
BCF-080	Bladder Filling Catheter 8fg	25
BCF-010	Bladder Filling Catheter 10fg	25
BCF-012	Bladder Filling Catheter 12fg	25
DLBC-060	Double Lumen Cystometry Catheter 6fg	25
DLBC-080	Double Lumen Cystometry Catheter 8fg	25
DLBC-010	Double Lumen Cystometry Catheter 10fg	25
RBC-045	Rectal Balloon Catheter 4.5fg	25
RBC-1045	Rectal Balloon Catheter 4.5fg + Introducer	25
RPC-145	Rectal Pressure Catheter 4.5fg	25
UPPD-060	Double Lumen UPP Catheter + Yconn.6fg	25
UPPS-080	Single Lumen UPP Catheter 8fg	25
CP006	Cystometry Set (6fg Double Lumen & Rectal Balloon Catheter)	25
CP007	Cystometry Set (8fg Double Lumen & Rectal Balloon Catheter)	25