

Grip-Lok[®] 3400L Large Universal Securement

Now available on FP10 Drug Tariff

Grip-Lok[®] 3400L is a versatile stabilisation and securement device for locking catheters, tubes and lines securely in place.

Hypoallergenic, breathable and latex-free skin contact adhesive

Provides superior securement for both vertical and horizontal lifting accidental line pulls

Universally works with tubes, lines or catheters sized 6 - 13mm in diameter

Soft and flexible fabric design improves patient comfort

Easy to apply, open, inspect and adjust

Sterile packed to allow application of the device close to the insertion site

Low profile with no pins or hard plastic parts

Can easily be placed and adjusted with gloved hands

Grip-Lok[®] 3400L secures chest tubes, Foley catheters or any tube, line or catheter sized between 6mm and 13mm in diameter.



Reference	NPC	Description	Tubing Range	Device Length	Units/Box
3400L	ELW298	Grip-Lok Large Universal Securement	18 - 40 Fr	15.9cm	100



BLUE BOX
MEDICAL LTD

Unit 29, New Forest Enterprise Centre,
Chapel Lane, Totton, Southampton SO40 9LA



Grip-Lok is a trademark of TIDI Products, LLC



Tel: 023 8066 9000 Fax: 023 8066 9002
info@blueboxmedical.co.uk
www.blueboxmedical.co.uk



PS041

05/16 - B

FIXATION