

Auto PC[®] Automatic Pressure Infuser

PLEASE SEE OVERLEAF FOR SPECIFICATION

Wider range of pressure set-up:
100mmHg ~ 400mmHg

Use of various reservoirs:
250ml, 500ml, 1,000ml etc.
Able to use normal fluid or blood bag

Safety alarm:
over pressure
(+20mmHg),
low pressure
(-20mmHg),
sensor out,
air leakage

Noise
reduction



Single type



Rechargeable type



Double type

Various applications:

- Constant stable pressure for BP monitoring
- Fast infusion for blood transfusion

No need for extra power pack

Max pressure reaching time: 20 sec


Can be used alone as a
pressure infuser and also
together with AutoMer II

BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002
info@blueboxmedical.co.uk
www.blueboxmedical.co.uk

ACEMEDICAL

PATIENT CARE



PS059

03/17 - A

AUTO PC[®] AUTOMATIC PRESSURE INFUSER

SPECIFICATION

Classification	Class II b
Power Source	DC 12V, 2A
Pressure Range	100 ~ 400mmHg (Steps: 5mmHg)
Alarm	Over Pressure (+20mmHg), Low Pressure (-20mmHg) Pressure Sensor Out, Air Leakage
Weight	2.2Kg (Single), 4.9Kg (Double)
Dimensions	Single Type: 190(W) x 385(H) x 100(D) mm Double Type: 458(W) x 385(H) x 105(D) mm
Conformity	EN60601-1, EN60601-1-2

BLUE BOX
MEDICAL LTD 

Tel: 023 8066 9000 Fax: 023 8066 9002

info@blueboxmedical.co.uk

www@blueboxmedical.co.uk

ACEMEDICAL